



## 1 EPOXY PEGBOARD

### sizes

- 24" wide x 30" high - 39 pegs - EPB2430
- 30" wide x 30" high - 46 pegs - EPB3030
- 30" wide x 36" high - 56 pegs - EPB3036

- material: 1" thick black (or gray) epoxy resin with beveled edges
- installation: 4 screw holes provided for wall attachment
- drip trough: 2" d stainless steel with 36" clear flexible drain tube
- some assembly required

## 2 MIXING FAUCET

### mixing faucet (with VB)

L414-8VB-55-BO140S

- mixing faucet: deck-mounted, with vacuum breaker for hot and cold water
- handles: forged brass four-arm handles with color-coded index discs
- finish: forged brass with polished chrome plated finish

### faucet

valve units: self-contained compression valve with replaceable stainless steel seats

gooseneck: 8" rigid/swing convertible gooseneck

inlet: 1-3/16" male shank with (2) 3/8" OD flexible copper tubes; furnished complete with locknut and washer

outlet: 3/8" NPS female outlet; included ten serration hose end or aerator

vacuum breaker (VB): integral vacuum breaker to prevent backflow; breaker has replaceable seat and ultra-light float cup to resist spilling at low flow

## 3 BLACK EPOXY RESIN

### work surfaces

- 30" x 36" - EPX3036
- 30" x 48" - EPX3048
- 30" x 60" - EPX3060
- 30" x 72" - EPX3072
- 30" x 96" - EPX3096
- 30" x 100" - EPX30100
- 36" x 100" - EPX36100

### sink top\*

- 30" x 50" - EPX3050E
- 30" x 60" - EPX3060E
- 30" x 72" - EPX3072E

### backsplash

- 4" x 60" - EPX0460
- 4" x 100" - EPX04100

### sidesplash

- 4" x 29" - EPX0429

### work surface

1" thick monolithic epoxy resin

edges: 1/8" bevel with drip groove on underside

splash: applied

adhesive: two-part epoxy provided as needed

other sizes can be readily fabricated; please consult your representative for lead times

\*sinktop E: E indicates that the top includes a black D55 dop-in epoxy sink (25"x15"x10"), cutout and outlet; sinks are centred right to left in the work surface

## GRAY EPOXY RESIN

**work surfaces**

30" x 36" - EPX3036  
 30" x 48" - EPX3048  
 30" x 60" - EPX3060  
 30" x 72" - EPX3072  
 30" x 96" - EPX3096

**sink top\***

30" x 50" - EPX3050E

**backsplash**

4" x 60" - EPX0460  
 4" x 100" - EPX04100

**sidesplash**

4" x 29" - EPX0429

- 1" thick monolithic epoxy resin
- edges: 1/8" bevel with drip groove on underside
- adhesive: two-part epoxy provided as needed

**work surfaces**  
 splash: applied

other sizes can be readily fabricated;  
 please consult your representative for  
 lead times

**E - epoxy sinks**  
**S - stainless stinks**

\*sink tops

sinktop E: E indicates that the top includes a black d55 dop-in epoxy sink (25"x15"x10"), cutout and outlet; sinks are centred right to left in the work surface

sinktop S: S indicates SLXD2225-A-GR-JX stainless steel drop-in sink (22"x16"x12"), cutout and strainer; sink ledge can accommodate a quick-ship mixing faucet and eyewash as well as single faucet; left 2 holes come with caps; sinks are centred right to left in the work surface

see page 20 for black epoxy resin work surfaces

## BLACK PHENOLIC RESIN

**work surfaces**

30" x 36" - PHN3036  
 30" x 48" - PHN3048  
 30" x 60" - PHN3060  
 30" x 72" - PHN3072  
 30" x 96" - PHN3096  
 30" x 120" - PHN30120

**sink top\***

30" x 60" - PHN3060E  
 30" x 72" - PHN3072E  
 30" x 60" - PHN3060S

**backsplash**

4" x 120" - PHN04120

**sidesplash**

4" x 29" - PHN0429

- 1" thick laminated kraft paper finished on all 4 edges and both faces
- edge: 1/8" bevel

## GRAY PHENOLIC RESIN

**work surfaces**

30" x 36" - PHN3036  
 30" x 48" - PHN3048  
 30" x 60" - PHN3060  
 30" x 72" - PHN3072  
 30" x 96" - PHN3096  
 30" x 120" - PHN30120

**sink top\***

30" x 60" - PHN3060E  
 30" x 72" - PHN3072E  
 30" x 60" - PHN3060S

**backsplash**

4" x 120" - PHN04120

**sidesplash**

4" x 29" - PHN0429

- 1" thick laminated kraft paper finished on all 4 edges and both faces
  - edge: 1/8" bevel
  - please note: gray phenolic will have a black edge
- \*gray phenolic with gray core can be ordered for an additional fee

## BLACK PLASTIC LAMINATE

**work surfaces**

30" x 36" - PLM3036  
 30" x 48" - PLM3048  
 30" x 60" - PLM3060  
 30" x 72" - PLM3072  
 30" x 96" - PLM3096  
 30" x 120" - PLM30120

**backsplash**

4" x 120" - PLM04120

**sidesplash**

4" x 29" - PLM0429

- p-lam work surfaces: 1" thick medium density particleboard with 1/16" thick laminate on both faces finished with a durable 3mm PVC edge on all 4 edges
- edges: 1/8" radius
- splash: applied
- color: black (matching edgeband)
- standard grade

## GRAY PLASTIC LAMINATE

**work surfaces**

30" x 36" - PLM3036  
 30" x 48" - PLM3048  
 30" x 60" - PLM3060  
 30" x 72" - PLM3072  
 30" x 96" - PLM3096  
 30" x 120" - PLM30120

- p-lam work surfaces: 1" thick medium density particleboard with 1/16" thick laminate on both faces finished with a durable 3mm PVC edge on all 4 edges
- edges: 1/8" radius
- splash: applied
- color: gray (matching edgeband)
- standard grade

## STAINLESS PEGBOARDS



### sizes

18" wide x 24" high - 15 pegs - V1824  
 24" wide x 30" high - 32 pegs - V2430  
 30" wide x 36" high - 40 pegs - V3036

- material: 20 gauge, type 304 stainless steel, polished to #4 finish
- installation: wall mount bracket shipped with each pegboard
- 6" long 1/2" diameter white polypropylene pegs
- cleaning: safe to autoclave
- integrated drip tray

see page 20 for epoxy pegboards

## EPOXY SINK



### sink

D55

- drop-in
- 25" w x 15" h x 10" d
- cutout and outlet
- colors: black and gray



### outlet

OUTL-S03R-001

### sink strainer

STRN-01-001

- for use with D55 sink
- colors: black and gray

## STAINLESS STEEL SINK



### sink

SLXD-2225

- 22" w x 16" h x 12" d
- cutout and strainer

## ACRYLIC SPLASHGUARD



### sizes

30" x 12" x 0.5" - ASGN3012-SS29-SS8  
 30" x 18" x 0.75" - ASGN3018-SS29-SS14

- notched with stainless steel U-channel
- eased and flame-polished edges
- 2" radius upper corner

